

Added Comfort and Safety with our Fall Detection Help Button!

An additional safeguard for you, our Fall Sensor help button is programmed with our standard medical alert systems (cellular and landline options available) and was designed for situations when you fall and are unable to push the help button.

How it Works?

When a fall or accident occurs and no sudden movement is detected within 10 seconds, an "unconsciousness" alarm will be raised and sent to the Care Center to get you the help you need.

Direct Link will then activate the unit and talk to you to ensure that you are safe. In the event of an emergency, our Care Center will contact a caregiver, next door neighbor, or emergency personnel when appropriate, so that the most appropriate actions can be taken.

Fall Detection Button Features:

- Water resistant
- Battery lasts two years with normal usage
- Alerts are sent to the Care Center when a low battery is detected.
- Help button can be pressed from up to 600 feet from base unit



Used in conjunction with our standard medical alert systems (cellular and landline options available).

***Call to learn more!
(630) 323-7231***

